

# **EXPERIMENTS PROCEDURE**

## **Full Explanation on \*TITRATION\***

Titration is a process whereby one solution(acid) is added gradually to another solution(base) until end point is reached.

End point is the point where the acid neutralizes the base or viceversa.

During titration, a phenolphthalein indicator is used.

Phenolphthalein is pink in colour.

When added to a base,the mixture

turns pink but when phenolphthalein is added to an acid the mixture turns red.

During the titration experiment, one solution must have its concentration known. A solution with known concentration is called a \*Standard solution\*

Usually, (but not always), Sodium Hydroxide ( $\text{NaOH}$ ) has its concentration known.

And Hydrochloric acid ( $\text{HCl}$ ) has unknown concentration.

After titration we need to calculate

the molarity or concentration of the unknown that's how much has been used to neutralize the standard solution.

Now the whole experiment.

Suppose Sodium hydroxide ( $\text{NaOH}$ ) is our standard solution and Hydrochloric acid is our unknown concentration.

i) Firstly, we add eg 20ml of  $\text{NaOH}$  (base) into a conical flask. Record the volume (eg 20ml)

ii). Using a dropper, we add 2 or 3

drops of Phenolphthalein indicator into the conical flask with NaOH. The mixture will be pink in colour.

iii). Then we measure 20ml of HCl and we are told to add to the mixture of NaOH and Phenolphthalein indicator.

The adding should be done slowly by slowly or gradually while shaking the contents.

Stop adding the acid until the end point is reached thus the mixture has changed from pink into colourless.

This is the point where the acid has neutralized the base.

iv) Record the volume used to neutralize the base ( e.g. if it was 20ml and mwatsala 13ml then 7ml has bn used in the reaction)

v) Calculate the concentration of unknown using the formula

$C_1V_1=C_2V_2$ , where

C=concentration or molarity and V=Volume.

**\*Note\***

During titration you maybe asked to write the balanced equation for the

reaction so it is



This one is already balanced.

Some questions also follow after the experiment like

Source of error and how to reduce.

Error source could overshooting.

.....

**\*FOOD TESTS\***

Food substances required: solution

of glucose, sucrose and starch;  
ground-up.

suspension in water of  
pea,bean,carrot, and grape  
;milk,castor -oilbeans, and  
cookingoil;

pieces of potato, bread ,and boiled  
white of egg

( Note. Sometimes you are provided  
with the real food stuff koma Pena  
you are given the food solution but  
the procedure is the same )

1- \*Test for Starch\*  
-Method

place about 3cm of  
starch solution in a  
container

add 2-3drops of iodine  
solution

observe colour change

-Results

blue-black  
color

means( presence of  
starch

)

browncolor means ( absence of  
starch

)



## 2- \*Test for Protein\*

### - Method

place about  $3\text{cm}^3$  of the food solution in a test tube

carefully add 4-5 drops of biuret solution.

### -Results

purple color means( presence of protein

)

blue colour means (absence of protein

)

## ii- \*Test for Protein\*

### - Method

placea bout  $3\text{cm}^3$  of food  
solution in a test tube

add about  $3\text{cm}^3$  of sodium  
hydroxidesoln.and mix  
well

-Add 1or 2drops of  
coppersulphate solution mix well and  
observe  
any colour change

NB:if the colour does not change  
immediately let the solution stand  
for 15minutes and then observe

again

## -Results

violet color (presence of protein)

## 3- \*Test for Lipids\*

### -Method

place a drop of oil/food substance on a piece of paper (filter paper)

examine the paper while holding up against light

### -Results

translucent  
or grease  
spot  
means (presence of  
lipids  
)  
not  
transparent  
means ( absence of  
lipids)

## ii- \*Test for Lipids\*

2drops of cooking oil are  
thoroughly shaken with  
about 5cm<sup>3</sup> of ethanol in a

dry test tube until the fat dissolves.

the alcoholic solution is poured into a test tube containing a few water

(NB: this is an ethanol/emulsion test  
)

-Results

cloudy

white

emulsion

means (presence of fats/lipids

)

if the liquid

stay clear

means ( absence of  
fats

)

4- \*Test for Reducing  
sugars\*

(monosaccharides .eg  
glucose

and

maltose)

-Method

place about  $3\text{cm}^3$  of the

food solution in a test tube

add about 2 drops of

benedict solution

heat gently for few

minutes

(NB: the mouth of the test-tube

must be pointed away from

people as the solution tends to bump

out of the tube)

-Results

A yellow or

orange or

brick red

color

means( presence of  
reducing  
sugar  
)

blue colour means (absence of  
reducing  
sugar  
)

## ii- \*Test for Reducing Sugars\*

### - Method

place about  $3\text{cm}^3$  of the  
food solution in a test tube

Add about  $2\text{cm}^3$  of  
Fehling's A solution followed



by the same Fehling'sB  
solution.

heat gently for few  
minutes

-Results

Ayellow or  
orange or  
brickred  
color

means( presence of  
reducing  
sugar),,,

# \*EFFECTS OF PRACTICE ON HITTING THE TARGET\*

## \*Procedure\*

- 1-Put a mark on the wall
- 2-measure a distance of 10 m from the mark and draw a line
- 3- identify a person and give him 10 small objects such as darts or stones
- 4-let the person stand on the line and throw one object at a time towards the mark
- 5-count the number of objects that hit the mark and record it
- 6- repeat the same activity with the

same person for three more trials .

7-compare the number of objects that hit the mark in all trials .

**\*Expected results\***

More objects will hit the mark in the last trial than the first trial

**\*Conclusion\***

Practice increases chances of hitting the target

**\*EFFECTS OF TIME OF THE DAY  
ON MEMORISING LIST OF  
WORDS\***

## \*Procedure\*

- 1- write a list of different words of same length and level of difficult on three chart papers
- 2- Identify the person who will recite the words on the paper.
- 3-flash the first paper to a person for a minute in the morning and hide it
- 4-Ask the person to recite the words of the hidden paper.
- 5-Record the number of words that are recited correctly by the person
- 6- Repeat the same activity with the same person using the second and third paper at noon and in the afternoon respectively

7-Compare the results obtained in three trials .

**\*Expected results\***

More words will be corectly recited in the morning than at noon and in the afternoon

**\*Conclusion\***

People memorise more words in the morning than at noon and in the afternoon

**\*EXPERIMENT TO SHOW THAT A**

# GERMINATING BEAN SEED CONTAINS ENZYMES THAT DIGEST STARCH\*

## \*Procedure\*

1. crush the germinating bean seed in a mortar
2. add water to obtain(enzyme) extract
- 3.put starch soln. in test tube A and B
- 4.in a test tube A add extract ,but leave test tube B in intact
- 5.leave both test tube to stand for some time
- 6.later add drops of iodine soln. to

both test tubes and observe color change

### **\*Results\***

In test tube A-there will be brown color. This show that there's no starch present

In test tube B –you will observe blue black color. This change indicate presence of starch

### **\*conclusion\***

Germinating bean seeds contain an enzyme

which digested starch in test

All de best

@Mr President

Lighton D Malunga

0882838030